

ASLP

Compact Hand System

1.5

2.0

2.4



2.0 mm ASLP Compact Hand System

Indications :

- Fractures of the phalanges, metacarpals, distal radius and wrist bones.
- Osteotomies and arthrodeses of the interphalangeal joints.



2.0 Cortex Screw

Length 6 to 14 mm (1 mm Increment)
Length 16 to 34 mm (2 mm Increment)
Length 38 mm



2.0 Locking Screw

Length 6 to 14 mm (1 mm Increment)
Length 16 to 30 mm (2 mm Increment)



2.0 Adaption Plate



2.0 Mini H-Plate



ASLP - 2.0 Adaption Plate



ASLP - 2.0 Condylar plate



ASLP - 2.0 Rotation Correction Plate



ASLP - 2.0 Straight Plate




ASLP - 2.0 Y - Adaption Plate




ASLP - 2.0 T - Plate



ASLP - 2.0 T - Adaption Plate



2.0 mm System



2.4 mm System