

# ASLP-2.4/2.7 Distal Ulna Plate



## ASLP-2.4/2.7 Distal ulna plate

| PT362.03 - **3 Hole**

| PT362.04 - **4 Hole**

| PT362.05 - **5 Hole**

### Indication:

- Indicated for distal ulna metaphyseal fracture.

### Feature:

- Precontoured universal plate for distal ulna bone available with 4 nos plate head locking holes and 3,4 and 5 shaft combi holes variance.
- Distal ulna plate accommodates 2.4 mm Locking screw in the Plate head and 2.4 mm Locking, 2.4 & 2.7cortex screws in combi hole in the shaft region.
- Available in Titanium Alloy.

