



Headless Compression Screws




Indication:

- Fixation of small bone and the fragment of small bone.
- Intra articular fracture of tarsals, metatarsals, carpals, metacarpals and extra articular fractures.
- Fracture of the patella, ulna and radial styloid.
- For scaphoid bone fixation.

* Threads are of equal lead (i.e they advance the same distance for each revolution)



Headless Compression Screws



∅ 2.4mm Headless Compression Screws, short thread. (Titanium)

| ST273.10 - L 10mm
| ST273.12 - L 12mm
| ST273.14 - L 14mm
| ST273.16 - L 16mm
| ST273.18 - L 18mm
| ST273.20 - L 20mm
| ST273.22 - L 22mm
| ST273.24 - L 24mm
| ST273.26 - L 26mm
| ST273.28 - L 28mm
| ST273.30 - L 30mm
| ST273.32 - L 32mm
| ST273.34 - L 34mm
| ST273.36 - L 36mm
| ST273.38 - L 38mm
| ST273.40 - L 40mm

∅ 2.4mm Headless Compression Screws, long thread. (Titanium)

| ST274.18 - L 18mm
| ST274.20 - L 20mm
| ST274.22 - L 22mm
| ST274.24 - L 24mm
| ST274.26 - L 26mm
| ST274.28 - L 28mm
| ST274.30 - L 30mm
| ST274.32 - L 32mm
| ST274.34 - L 34mm
| ST274.36 - L 36mm
| ST274.38 - L 38mm
| ST274.40 - L 40mm



∅ 3.00 mm Headless Compression Screws.

| ST105.12 - L 12mm | ST105.28 - L 28mm
| ST105.14 - L 14mm | ST105.30 - L 30mm
| ST105.16 - L 16mm | ST105.32 - L 32mm
| ST105.18 - L 18mm | ST105.34 - L 34mm
| ST105.20 - L 20mm | ST105.36 - L 36mm
| ST105.22 - L 22mm | ST105.38 - L 38mm
| ST105.24 - L 24mm | ST105.40 - L 40mm
| ST105.26 - L 26mm

