




ACP™ - Anterior Cervical Plating System




Sizes :

Material	Level	Length (mm)
Titanium Alloy	I	19, 21, 23, 25, 27.5, 30, 32.5
	II	35, 37.5, 40, 42.5, 45, 47.5, 50
	III	52.5, 55, 57.5, 60, 62.5, 65, 67.5, 70, 72.5



Adequate window for graft visibility.

Integrated Locking Mechanism.



- Integrated locking mechanism enabled with ¼ twist of a driver prevents screw backing out.



Self- Drilling, Self Tapping Screws

- Reduces the need for intra-operative surgical steps.

Sizes :

Description	Height (mm)
4.0 mm Cervical Screw	11, 13, 15, 17, 19
4.5 mm Cervical Revision Screw	11, 13, 15, 17, 19



ACP™ - Anterior Cervical Plating System



Basic Instrument Set

Sr. No.	Description	Code
1	Tap	01319-N
2	Drill bit 2.5 mm caliberated stopper	01320-N
3	Drill bit 2.5 mm caliberated stopper	01321-N
4	Drill bit 2.5 mm	01322-N
5	Hexagonal driver shaft	01323-N
6	Hexagonal driver shaft	01324-N
7	Cervical plate holding forcep	01325-N
8	Locking spanner for secure lock	01326-N
9	Screw holding sleeves	01327-N
10	Curette	01328-N
11	Plate bender	01329-N
12	Pin distractor	01330-N
13	Starter Awl with adjusting stopper	01331-N
14	Temporary pin - Short	01332-N
15	Temporary pin - Long	01333-N
16	Staight AO QC handle	01334-N
17	Drill guide for 2.5 mm drill bit	01335-N
18	Temporary pin holder for short pins	01336-N
19	Temporary pin holder for Long pins	01337-N
20	Plate box	-

INSTRUMENT SET

