


ALC™ -Plif Cage




Features:

- Tantalum marker for implant visualization.
- Large area for bone graft.
- Convex shape to match natural anatomical design.
- Bullet nose for easy insertion.


Sizes :

Material	Degree	Width	Height (mm)
PEEK	0°	22 / 26	7, 8, 9, 10, 11, 12, 13, 14
Titanium	0°	22 / 26	7, 8, 9, 10, 11, 12, 13, 14



- Toothed upper & lower surfaces for better anchoring.

ALC™ -Plif Cage



Instrument Set

