

# ASLP 3.5-5.0 Distal Lateral Condylar Plate



## Indication:

1. Designed for supplementary fixation in femoral lateral condyle comminuted fractures, in combination with a Distal Femur Nail.

## Feature:

1. 7-hole locking plate system designed to integrate with the Distal Femur Nail, providing enhanced angular stability.
2. Anatomically precontoured plate for precise fit and alignment.
3. Improves implant anchoring in cases of poor bone quality and facilitates early patient mobilization.
4. Two 5.0mm locking screw/5.0mm advanced IM nail screw ensures stable plate-nail combination.
5. Additional peripheral support with five 3.5mm locking screws.


# ASLP 3.5-5.0 Distal Lateral Condylar Plate



**ASLP - 3.5/5.0 Distal Lateral Condylar Plate**

PT480.07.01  
Right | 7 Holes

PT480.07.02  
Left | 7 Holes



### 3.5mm Locking Screw, Self Tapping

| Length 8 mm to 50 mm (2 mm increment)  
Also available in *ST032.08 to ST032.50*

| Length 55 mm to 95 mm (5 mm increment)  
Also available in *ST032.55 to ST032.95*


| Length 35 mm & 45 mm  
Also available in *ST032.35 & ST032.45*



### 5.00mm Locking Screw, Self Tapping

| Length 18 mm to 50 mm (2 mm increment)  
Also available in *ST033.18 to ST033.50*

| Length 55 mm to 95 mm (5 mm increment)  
Also available in *ST033.55 to ST033.95*



### 5.00mm Advanced IM Nail Screw

| Length 26 mm to 70 mm (2 mm increment)  
Also available in *ST461.26 to ST461.70*

| Length 75 mm to 100 mm (5 mm increment)  
Also available in *ST461.75 to ST032.100*

