

ASLP™ - 3.5, Clavicle Hook Plates.

Elongated combihole facilitates adjustment of plate on bone.


Rounded edge and posterior hook offset.

Anatomical preconture plates with left and right variance.

Low profile distal head prevents soft tissue impingement

12° anatomic prebend enables accurate placement.

Various hook depth
12 mm, 15 mm and 18 mm.



ASLP™ - 3.5, Clavicle Hook Plates.

ASLP - 3.5, Clavicle Hook Plates

Titanium	Shaft holes		Hook depth
PT120.12.04.1.12	4	right	12
PT120.12.04.2.12	4	left	12
PT120.12.04.1.15	4	right	15
PT120.12.04.2.15	4	left	15
PT120.12.04.1.18	4	right	18
PT120.12.04.2.18	4	left	18
PT120.12.05.1.12	5	right	12
PT120.12.05.2.12	5	left	12
PT120.12.05.1.15	5	right	15
PT120.12.05.2.15	5	left	15
PT120.12.05.1.18	5	right	18
PT120.12.05.2.18	5	left	18
PT120.12.06.1.12	6	right	12
PT120.12.06.2.12	6	left	12
PT120.12.06.1.15	6	right	15
PT120.12.06.2.15	6	left	15
PT120.12.06.1.18	6	right	18
PT120.12.06.2.18	6	left	18
PT120.12.07.1.12	7	right	12
PT120.12.07.2.12	7	left	12
PT120.12.07.1.15	7	right	15
PT120.12.07.2.15	7	left	15
PT120.12.07.1.18	7	right	18
PT120.12.07.2.18	7	left	18
PT120.12.08.1.12	8	right	12
PT120.12.08.2.12	8	left	12
PT120.12.08.1.15	8	right	15
PT120.12.08.2.15	8	left	15
PT120.12.08.1.18	8	right	18
PT120.12.08.2.18	8	left	18

