

# ASLP™-3.5


## T-Plate Oblique Angled

### Indication :

- Intended use for fixation of fractures of Distal Radius & Tibia bone.

### Features :

- Plates comprise of 3 locking holes in the head and combi holes in the shaft.
- Anatomical plates with right and left options.
- Accommodates 3.5 mm locking screws and 3.5 mm cortex screw.
- Available in Titanium alloy.



### 3.5 T-Plate Oblique Angled

Titanium	Head Holes	Shaft Holes	R/L
PT237.33.01	3	3	Right
PT237.34.01	3	4	Right
PT237.35.01	3	5	Right
PT237.33.02	3	3	Left
PT237.34.02	3	4	Left
PT237.35.02	3	5	Left