


# ASLP-2.4/2.7/3.5 Posterior Fibula Plate



## Indication:

- For fixation of Distal Fibula fracture.
- Intended use for Spiral fibula fracture.
- For treating Non unions and Mal union of fractures.

## Feature:

- Anatomical precontoured locking plate matches to anatomical geometry of fibula bone.
- Plates are available in left and Right version with 3,4,5,6,7 and 8-hole variance.
- Accommodate 2.4mm/2.7 mm Locking screw in plate head and 3.5 mm Locking and cortical screws in the plate shaft.
- Elongated holes in the shaft region facilitates plate adjustment on the bone.



### ASLP-2.4/2.7/3.5 Posterior Fibula Plate

PT371.03.01 Right   3 Hole	PT371.03.02 Left   3 Hole
PT371.04.01 Right   4 Hole	PT371.04.02 Left   4 Hole
PT371.05.01 Right   5 Hole	PT371.05.02 Left   5 Hole
PT371.06.01 Right   6 Hole	PT371.06.02 Left   6 Hole
PT371.07.01 Right   7 Hole	PT371.07.02 Left   7 Hole
PT371.08.01 Right   8 Hole	PT371.08.02 Left   8 Hole

