


# 4.0 mm partially threaded Cannulated Cancellous Screw



## Indication:

- Fixation of tarsal and metatarsal fractures.
- Fixation in osteotomies of metatarsal and phalangeal.
- Indicated in tarsometatarsal and metatarsophalangeal arthrodesis.
- For fixation of ligaments.
- Used for hallux valgus correction.

## Feature:

- 4.0 mm CC screws possess reverse cutting flute.
- Fixation can be done with and without washers.
- 1.5 mm Guide wire is used for 4.0 mm CC screws.

### Ø4.0 mm Cannulated Screws, partially threaded

ST415.20 - L 20mm	ST415.48 - L 48mm
ST415.22 - L 22mm	ST415.50 - L 50mm
ST415.24 - L 24mm	ST415.52 - L 52mm
ST415.26 - L 26mm	ST415.54 - L 54mm
ST415.28 - L 28mm	ST415.56 - L 56mm
ST415.30 - L 30mm	ST415.60 - L 60mm
ST415.32 - L 32mm	ST415.64 - L 64mm
ST415.34 - L 34mm	ST415.68 - L 68mm
ST415.36 - L 36mm	ST415.72 - L 72mm
ST415.38 - L 38mm	ST415.74 - L 74mm
ST415.40 - L 40mm	ST415.76 - L 76mm
ST415.42 - L 42mm	ST415.78 - L 78mm
ST415.44 - L 44mm	ST415.80 - L 80mm
ST415.46 - L 46mm	