

# CTL™ -Mesh Cage System



## Indication:

- Intended used as vertebral body replacement in Thoraco-Lumber spine to replace the damaged, collapsed vertebral body due to tumor or trauma injury.

## Features:

- Rounded edges prevent the soft tissue irritation.
- Enable easy insertion.
- Can be used as cut to desire specific size.

## Sizes :

Diameter (mm)	Length (mm)										
	15 mm	20 mm	25 mm	30 mm	35 mm	40 mm	45 mm	50 mm	60 mm	70 mm	100 mm
10	•	•	•	•				•	•	•	•
12	•	•	•	•	•	•					•
14		•	•	•	•	•	•				
16			•	•	•	•	•				

