


MS-2 Fenestrated Poly Axial Screw




Indication:

- Degenerative Disc disease
- Spondylolisthesis
- Tumor and Trauma
- Spinal stenosis
- Non-fusion in previous surgery.
- Pseudo arthrosis
- Indicative for highly compromised bone quality

Feature and benefit:

- Cannulation of the screw is 1.8 mm
- Total 6 fenestrations in cannulated pedicle screw



Size chart and reference product code

Screw	Screw Length (mm)		
	35 mm	40 mm	45 mm
Diameter			
5.5 mm	•	•	•
6.2 mm	•	•	•
7.0 mm	•	•	•

