



MIS – II Pedicle Screw System



MIS – II Pedicle Screw System



Indication:

- For the use in Lumber fusion
- Spinal deformities
- For the Thoraco-Lumber Degenerative Pathologies
- Tumors and Trauma
- Spondylolisthesis

Features:

- Integrated Extended screw tab provide secured pathway of the rod implantation.
- Surgeons inspired instrumentation provides intraoperative assurance and reliability.
- Useful in multi-level fixation.
- Star head inner nut recess adequate tightening and reduces wear.
- Ergonomic instrument design provides ease of doing percutaneous surgery.

Benefits of MIS spine surgery :



- Decreased Pain
- Minimal Blood loss
- Faster Recovery
- Shorter hospital stay
- Less surgical complication

Color coding implants:



MIS-II Poly axial pedicle screw

ø5.5 Length 35 mm <i>Available in ST303.55.35</i>	ø6.5 Length 40 mm <i>Available in ST303.65.40</i>
ø5.5 Length 40 mm <i>Available in ST303.55.40</i>	ø6.5 Length 45 mm <i>Available in ST303.65.45</i>
ø5.5 Length 45 mm <i>Available in ST303.55.45</i>	ø6.5 Length 50 mm <i>Available in ST303.65.50</i>
ø6.0 Length 35 mm <i>Available in ST303.60.35</i>	ø7.0 Length 35 mm <i>Available in ST303.70.35</i>
ø6.0 Length 40 mm <i>Available in ST303.60.40</i>	ø7.0 Length 40 mm <i>Available in ST303.70.40</i>
ø6.0 Length 45 mm <i>Available in ST303.60.45</i>	ø7.0 Length 45 mm <i>Available in ST303.70.45</i>
ø6.0 Length 50 mm <i>Available in ST303.60.50</i>	ø7.0 Length 50 mm <i>Available in ST303.70.50</i>
ø6.5 Length 35 mm <i>Available in ST303.65.35</i>	



MIS-II Inner

IT305



MIS-II Contoured rod




ø5.5 Length 30 mm <i>Available in RT304.55.30</i>	ø5.5 Length 75 mm <i>Available in RT304.55.75</i>
ø5.5 Length 35 mm <i>Available in RT304.55.35</i>	ø5.5 Length 80 mm <i>Available in RT304.55.80</i>
ø5.5 Length 40 mm <i>Available in RT304.55.40</i>	ø5.5 Length 85 mm <i>Available in RT304.55.85</i>
ø5.5 Length 45 mm <i>Available in RT304.55.45</i>	ø5.5 Length 90 mm <i>Available in RT304.55.90</i>
ø5.5 Length 50 mm <i>Available in RT304.55.50</i>	ø5.5 Length 95 mm <i>Available in RT304.55.95</i>
ø5.5 Length 55 mm <i>Available in RT304.55.55</i>	ø5.5 Length 100 mm <i>Available in RT304.55.100</i>
ø5.5 Length 60 mm <i>Available in RT304.55.60</i>	ø5.5 Length 105 mm <i>Available in RT304.55.105</i>
ø5.5 Length 65 mm <i>Available in RT304.55.65</i>	ø5.5 Length 110 mm <i>Available in RT304.55.110</i>
ø5.5 Length 70 mm <i>Available in RT304.55.70</i>	ø5.5 Length 115 mm <i>Available in RT304.55.115</i>



MIS-II Straight rod

ø5.5 Length 90 mm <i>Available in RT359.55.90</i>	ø5.5 Length 140 mm <i>Available in RT359.55.140</i>
ø5.5 Length 95 mm <i>Available in RT359.55.95</i>	ø5.5 Length 145 mm <i>Available in RT359.55.145</i>
ø5.5 Length 100 mm <i>Available in RT359.55.100</i>	ø5.5 Length 150 mm <i>Available in RT359.55.150</i>
ø5.5 Length 105 mm <i>Available in RT359.55.105</i>	ø5.5 Length 155 mm <i>Available in RT359.55.155</i>
ø5.5 Length 110 mm <i>Available in RT359.55.110</i>	ø5.5 Length 160 mm <i>Available in RT359.55.160</i>
ø5.5 Length 115 mm <i>Available in RT359.55.115</i>	ø5.5 Length 165 mm <i>Available in RT359.55.165</i>
ø5.5 Length 120 mm <i>Available in RT359.55.120</i>	ø5.5 Length 170 mm <i>Available in RT359.55.170</i>
ø5.5 Length 125 mm <i>Available in RT359.55.125</i>	ø5.5 Length 175 mm <i>Available in RT359.55.175</i>
ø5.5 Length 130 mm <i>Available in RT359.55.130</i>	ø5.5 Length 180 mm <i>Available in RT359.55.180</i>
ø5.5 Length 135 mm <i>Available in RT359.55.135</i>	

Instrument Set



Miraclus Orthotech Pvt. Ltd.

Plant -I :

Plot No. 1112, Phase-3, G.I.D.C., Vatva, Ahmedabad - 382 445.
Phone : +91 79 2583 1332, 2583 1333
Email : contact@miraclus.com | www.miraclus.com

Plant -II :

Survey No. 106, Village Vavdi, Opp. Safex Industries,
Behind Elite Boarding School, Vavdi,
Ta. Kheda-387550, Gujarat, India

