



MS-II™ - Catalogue




Dual Lead, Top Loading
Pedicle Screw System.

MS-II spinal fixation system



52° Angulation enables fixation in asymmetrical pedicle geometry



Rod Capture Mechanism is designed to prevent rod migration upon initial insertion.



Ø5.7 mm Rod with Laser marked contour lines.

Poly Axial Screw Colour-coding



Ø 4.5



Ø 5.5



Ø 6.2



Ø 7.0


- Colour coding implants for easy identification.

MS-II spinal fixation system




Sizes :

Screw	Screw Length (mm)						
	25 mm	30 mm	35 mm	40 mm	45 mm	50 mm	55 mm
4.5 mm	•	•	•	•	•		
5.5 mm			•	•	•	•	•
6.2 mm			•	•	•	•	•
7 mm			•	•	•	•	•

- Dual Lead enables secure Anchorage in cortical and cancellous bone.
- Dual thread pattern ensures pull-out strength.
- Threaded blunt tip enables prompt bone-grip wit safety.



Instrument Set



Miraclus Orthotech Pvt. Ltd.

Plant address :

Plot No. 1112, Phase-3, G.I.D.C., Vatva, Ahmedabad - 382 445.
Phone : +91 79 2583 1332, 2583 1333
Email : contact@miraclus.com | www.miraclus.com

Plant - II address:

Survey no 106, village vavdi, opp safex industries
Bh elite boarding school, vavdi, ta kheda-387550