

Humerus Nail System

Indication:

- Proximal humerus fracture for two and three part.
- Humerus shaft fracture.
- Entry point: On the greater tubercle to avoid the damage of the axillary nerve.

Feature:


- Proximal bend of Humerus nail is 7°.
- Proximal Locking accommodates 5 mm Cancellous locking bolt in oblique direction for rotational Stability.
- Distal locking accommodates 3.5 Locking bolt.



Proximal Locking



Distal Locking



Available Humerus Nail Sizes:

Dia (mm)	Types	Length (mm)							
		160	180	200	210	220	240	260	280
6.5	Solid	160	180	200	210	220	240	260	280
7	Cannulated	160	180	200	210	220	240	260	280
7.5	Cannulated	160	180	200	210	220	240	260	280
8	Cannulated	160	180	200	210	220	240	260	280

Humerus Nail Screw and End cap:



Ø 3.5 mm Humerus Nail Locking Screw

| Length(mm) 18 to 50
(2 mm increment)



Ø 5.0 mm Humerus Nail Locking Screw


| Length(mm) 18 to 50
(2 mm increment)



Humerus nail end cap

| Length(mm) 0 to 15
(5 mm increment)

Instrument Set:



Miraclus Orthotech Pvt. Ltd.

Plant address :

Plot No. 1112, Phase-3, G.I.D.C., Vatva, Ahmedabad - 382 445.
Phone : +91 79 2583 1332, 2583 1333
Email : contact@miraclus.com | www.miraclus.com

Plant - II address:

Survey no 106, opp Safex Industries Bh
Elite Boarding School, Vavdi, Ta Kheda-387550