

PEDI PLATE



Indication:

- Correction of acquired deformities and congenital deformities in a pediatric patients.

Features:

- It can be very well placed with a minimally invasive approach.
- Early weight bearing is possible under surgeon's observation.
- Redirect the bone growth in long bones and gradually correct the angular deformities.
- Screw hole design allows up to 60° screw angulation.
- Low profile plate accommodates \varnothing 4.0 mm cannulated screws.



\varnothing 4.0 mm cannulated Screws

- | | |
|-------------------|-------------------|
| ST287.16 - L 16mm | ST287.28 - L 28mm |
| ST287.18 - L 18mm | ST287.30 - L 30mm |
| ST287.20 - L 20mm | ST287.32 - L 32mm |
| ST287.22 - L 22mm | ST287.34 - L 34mm |
| ST287.24 - L 24mm | ST287.36 - L 36mm |
| ST287.26 - L 26mm | |



Pedi Plate

- | |
|-------------------|
| PT288.12 - L 12mm |
| PT288.16 - L 16mm |