


ASLP™ - 3.5

Periarticular Proximal Humerus Plate



- Plate can be fixed 10-15 mm distally from Rotator cuff to prevent acromial impingement.
- Low profile plate enables fixation without impinging soft tissue.
- Plate can be used in Osteopenic bone condition.



ASLP™ - 3.5

Periarticular Proximal Humerus Plate



ASLP - 3.5, Periarticular Proximal Humerus Plate

Titanium	Shaft holes	
PT068.02.01	2	Right
PT068.02.02	2	Left
PT068.03.01	3	Right
PT068.03.02	3	Left
PT068.04.01	4	Right
PT068.04.02	4	Left
PT068.05.01	5	Right
PT068.05.02	5	Left
PT068.06.01	6	Right
PT068.06.02	6	Left