

ASLP™-4.5/5.0/6.5 Proximal Femur Plate



ASLP-4.5/5.0/6.5 Proximal femur plate

PT224.06.01 | 6 Holes Right PT224.06.02 | 6 Holes Left

PT224.08.01 | 8 Holes Right PT224.08.02 | 8 Holes Left

PT224.10.01 | 10 Holes Right PT224.10.02 | 10 Holes Left

PT224.12.01 | 12 Holes Right PT224.12.02 | 12 Holes Left

PT224.14.01 | 14 Holes Right PT224.14.02 | 14 Holes Left



Indication:

- Fracture of Trochanteric region:
 - | Multifragmentary pertrochanteric fracture.
 - | Intertrochanteric fracture.
 - | Trochanteric diaphyseal and transverse fracture.
- Ipsilateral shaft fracture.
- Metastatic fracture of proximal femur.
- Osteotomies of proximal femur.
- Fracture Fixation of osteopenic bone, non-unions and mal-unions.
- Fracture of medial cortex of femur.
- Subcapital femoral neck fracture.
- Subtrochanteric femur fracture.

Features:

- Combi -hole plate accommodate 4.5 cortical screws and 5.0 mm locking screws.
- Pre-contoured shaped anatomical plate.
- Proximal femur bone-plate fixation is done with size 6.5 mm cannulated locking screws with desire angulation.
- Narrow tip of the plate enable easy insertion.



4.5 mm Cortex screw, Self-Tapping

| Length 18 to 70 mm (2 mm increment)
ST029.18 to ST029.70

| Length 55 & 65 mm
ST029.55 & ST029.65



5.0 mm Locking screw, Self-Tapping

| Length 18 to 50 mm (2 mm increment)
ST033.18 to ST033.50


| Length 55 to 95 (5 mm increment)
ST033.55 to ST033.95



6.5 mm Cannulated Locking screw

| Length 30 to 120 mm (5mm Increment)
ST234.30 to ST234.120

ASLP™-4.5/5.0 Curved Broad Plate



Indication:

- ASLP - 4.5 / 5.0 mm curved broad plate is used for various fracture of femur i.e. periprosthetic fracture, and fracture fixation in osteopenic bone structure as well.
- Used for fixation of nonunions and malunions of fractured bones.

Features:

- Matches anatomical bow shape anterior-posteriorly of Femur bone.
- Limited contact between plate and bone interface enable blood vascularization.
- Accommodates 5 mm conventional self- tapping locking screw and special unicortical self- tapping blunt locking screws.

ASLP - 4.5 / 5.0 mm Curved Broad Plates

PT178.16
| 16 Holes

PT178.20
| 20 Holes

PT178.18
| 18 Holes

PT178.22
| 22 Holes

5.0 mm Periprosthetic locking screw

ST223.08
| Length 8 mm

ST223.14
| Length 14 mm

ST223.10
| Length 10 mm

ST223.16
| Length 16 mm

ST223.12
| Length 12 mm

ST223.18
| Length 18 mm

