

ASLP-2.7/3.5 Proximal Humerus GT Plate



ASLP-2.7/3.5 Proximal Humerus GT Plate

| PT360.05.01 - 5 Hole Right

| PT360.05.02 - 5 Hole Left

Indication:

- For the avulsion fracture of Greater Tuberosity of proximal humerus bone.

Feature:

- Low profile pre-contoured plate with 1.6 mm thickness, which also prevent soft tissue irritation.
- Without limiting shoulder motion, avoid the postoperative subacromial impingement on the articular surface.
- Plate head accommodate 3 numbers of 2.7mm locking screw which create an angular stable construct to enhance high pullout strength
- Shaft hole of the plate accommodate 3.5 locking screw & 3.5 cortex screw.
- Possesses 5 nos K-wire holes for temporary fixation of plate and reduction, which can be also utilized for suturing the ligaments and tendons.