

ASLP™-3.5 Proximal Humerus Plate



Indication:

- Intended used to fixation of fragmented fractures of the proximal humerus bone.
- Pseudarthrosis in the humerus bone.
- Osteotomies in the proximal humerus bone.

Features:

- There is provision of 9 proximal holes which accommodates 3.5 locking head screws to provide stability and grip the multi-fragment fractures even in osteoporotic bones.
- Small proximal holes for suturing to maintain fracture reduction.
- Undercuts on the plate head enables suturing conveniently.
- Elongated hole of the plate shaft area for easy plate adaptability.



ASLP-3.5, Proximal Humerus Plates

PT070.03	PT070.05
3 Holes	5 Holes

ASLP-3.5, Long Proximal Humerus Plates

PT071.05	PT071.08	PT071.12
5 Holes	8 Holes	12 Holes

PT071.06	PT071.10
6 Holes	10 Holes

