


ASLP™-4.5/5.0 Proximal Lateral Tibia Plate



Indication:

- Indicated in proximal shaft fractures.
- For the fixation of metaphyseal fractures.
- Treating intra-articular fractures.
- Also for fixation of periprosthetic fractures.

Features:

- Anatomical plate available in left and Right options.
- Various length and holes are available.
- Accommodates 4.5 mm cortical and 5.0 mm locking screws.

ASLP - 4.5/5.0 Proximal lateral tibia plate

PT096.03.01 Right, 3 Holes	PT096.03.02 Left, 3 Holes
PT096.04.01 Right, 4 Holes	PT096.04.02 Left, 4 Holes
PT096.05.01 Right, 5 Holes	PT096.05.02 Left, 5 Holes
PT096.07.01 Right, 7 Holes	PT096.07.02 Left, 7 Holes
PT096.09.01 Right, 9 Holes	PT096.09.02 Left, 9 Holes
PT096.11.01 Right, 11 Holes	PT096.11.02 Left, 11 Holes
PT096.13.01 Right, 13 Holes	PT096.13.02 Left, 13 Holes